**EnVision™ Detection Kit,  
Peroxidase/DAB, Rabbit/Mouse**

FOR RESEARCH USE ONLY.

**Code No:GK500705**

**Storage: Store vial at 2-8°C.**

***Intended Use***

ChemMate™ EnVision™ Detection Kit, Peroxidase/DAB, Rabbit/Mouse, is intended for use in immunocytochemistry together with the TechMate™ and Autostainer Instruments. The kit is employed in a two-step procedure where the first step is incubation of the tissue with an optimally diluted primary rabbit or mouse antibody, and the second step is incubation with the ChemMate™ EnVision™/HRP, Rabbit/Mouse (ENV) reagent of the kit. This reagent is a peroxidase-conjugated polymer, which also carries antibodies to rabbit and mouse immunoglobulins. The reaction is visualized by the ChemMate™ DAB+ Chromogen also included in the kit.

***Reagents Materials provided***

**Bottle A ChemMate™ DAKO EnVision™/HRP, Rabbit/Mouse (ENV)**

5 mL, ready-to-use

Dextran coupled with peroxidase molecules and goat secondary antibody molecules against rabbit and mouse immunoglobulins. In buffered solution containing stabilizing protein and preservative.

**Bottle B ChemMate™ Substrate Buffer**

2\*6.25 mL

Buffered solution containing hydrogen peroxide and preservative.

**Bottle C ChemMate™ DAB+ Chromogen**

1\*0.25mL, 50 x concentrated

3,3'-diaminobenzidine tetrahydrochloride in organic solvent. The colour of this reagent may vary from strong violet to colourless without having any influence on the performance of the kit.

***Precautions :***

1. **Bottle A**, ChemMate™ EnVision™/HRP, Rabbit/Mouse (ENV), contains material of animal origin and it cannot be excluded that trace amounts of human material could be present due to manufacturing procedures. As with any product derived from biological sources, proper handling should be used.
2. Do not expose **Bottle C**, ChemMate™ DAB+ Chromogen or the Substrate Working Solution (CHROM) to strong light.
3. **Bottle C**, ChemMate™ DAB+ Chromogen contains 1.5% 3,3'-diaminobenzidine tetrahydrochloride, and is labelled: Harmful.

**Storage :** Store the kit at 2-8 °C. The prepared Substrate Working Solution (CHROM) should be stored at 2-8 °C and used within 5 days.

**Quality Control :** Each staining run should include a known positive control specimen to ascertain a proper performance of all the applied reagents. If the positive control specimen fails to demonstrate positive staining, labelling of test specimens should be considered invalid. A negative control reagent should be used with each specimen to identify any non-specific staining. If non-specific staining cannot be clearly differentiated from the specific staining, the labelling of the test specimen should be considered invalid.

***Please contact us for more information about this product.***

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